PNEUMOTM
PEAK get your
own summit

MANUAL DE USO PNEUMO PEAK

MANUAL DE USO PNEUMO PEAK

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Modelo: OLV-10SH

Marca: PNEUMO PEAK

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Empresa: PNEUMO PEAK SAS

ÍNDICE

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PNEUMO[®]
PEAK
get your
own summit

Introduction

General Description of the Equipment

Principles of Altitude Simulation

Equipment Components

Installation Requirements

Safety and Warnings

Preparation Before Use

Operating Instructions

Altitude Programs and Levels

Physiological Monitoring

Session Completion

Cleaning and Maintenance

Storage

Troubleshooting

Technical Specifications

Contraindications

Professional Use and Liability

Warranty and Technical Service

Contact

1. INTRODUCTION

PNEUMO[®]
PEAK
get your
own summit

This manual describes the correct and safe use of the PNEUMO PEAK altitude simulation equipment, designed to recreate normobaric hypoxia conditions through an air filtration system that produces deconcentrated and/or purified oxygen by means of pneumatic pumps, equivalent to different altitude levels above sea level.

The purpose of this equipment is to support programs for:

Sports training

Rehabilitation

Health and longevity optimization

Altitude acclimatization

✦ Important:

Read this manual completely before using the equipment.

2. GENERAL DESCRIPTION OF THE EQUIPMENT

The altitude simulation equipment is a system that reduces the fraction of inspired oxygen (FiO_2), simulating altitudes from sea level up to 9,000 meters above sea level.

It may be used during:

Rest

Physical training

Therapeutic protocols supervised by qualified professionals

3. PRINCIPLES OF ALTITUDE SIMULATION

Altitude simulation is based on normobaric hypoxia, which consists of reducing the oxygen concentration of the air without modifying atmospheric pressure.

Potential benefits include:

Cardiovascular system stimulation

Improved metabolic efficiency

Hematological adaptations

Physical performance optimization

4. EQUIPMENT COMPONENTS

The equipment consists of the following elements:

Hypoxia/Oxygen generator – 1 unit

Power cable – 1 unit

Modular display with electronic control – 1 unit

50-liter air reservoir with electronic oxygen sensor – 1 unit

Air delivery hose – 2 units

HEPA generator filter – 1 unit

Full-face mask – 1 unit

Mask facial support – 1 unit

Mask valve – 1 unit

User manual – 1 unit

Oxygen sensor accessory cable – 1 unit

Display or control interface – 1 unit

Physiological sensor – 1 unit

5. INSTALLATION REQUIREMENTS

Clean, dry indoor environment.

Environmental operating conditions:

The equipment must operate indoors at temperatures between +10 °C and +35 °C, with relative humidity between 15% and 80%, without condensation.

Storage conditions:

The equipment may be stored at temperatures between –10 °C and +50 °C, in a dry place protected from direct sunlight and heat sources.

Electrical connection: 110V, 60Hz (power adapter for 220V use)

Stable and level surface

Minimum free space around the equipment: 1 m²

6. SAFETY AND WARNINGS

 Important warnings:

Do not use without professional supervision in individuals with medical conditions

Do not use if the equipment shows visible damage

Stop use in case of dizziness, nausea, or discomfort

Do not modify the equipment

7. PREPARATION BEFORE USE

Before each session:

Check all connections

Confirm cleanliness of the mask or breathing system

Turn on the equipment

Select the program or altitude level

Place monitoring sensors

8. OPERATING INSTRUCTIONS

Properly place the mask or air delivery system

Adjust the desired altitude level

Start the session

Monitor physiological parameters throughout the session

9. ALTITUDE PROGRAMS AND LEVELS

Simulated Altitudes and Approximate FiO_2 Equivalence

Program	Simulated Altitude (m)	FiO_2 (%)
1	2500	15.1
2	3400	13.4
3	4300	12.3

4	5400	11.3
5	6400	10.6
6	6900	10.3
7	7600	9.8
8	7950	9.6
9	8600	9.1
10	9050	8.8

✦ Protocols must be defined by a trained professional.

10. PHYSIOLOGICAL MONITORING

During use, it is recommended to monitor:

Oxygen saturation (SpO₂)

Heart rate

Subjective user perception (Borg scale recommended)

Safety thresholds must be established by the responsible professional.

11. SESSION COMPLETION

Gradually reduce altitude level

End the session

Remove the mask or breathing system

Record session data

Verify return to resting physiological values

12. CLEANING AND MAINTENANCE

Daily cleaning of the mask and equipment is required.

The main unit must be cleaned with a soft cloth slightly dampened with water and neutral detergent. Do not use corrosive products, solvents, or excessive liquids.

The mask must be cleaned after each use. Wash with warm water and neutral detergent, followed by disinfection with 70% isopropyl alcohol or 3% hydrogen peroxide. Do not use chlorine, excessive heat, or aggressive sterilization methods. Allow the mask to dry completely before reuse.

Hoses must be cleaned with warm water and neutral detergent, thoroughly rinsed, and air-dried.

Disinfection according to center protocol

Periodic filter inspection

Technical maintenance every 6 months

13. STORAGE

Store in a dry place

Avoid direct sunlight exposure

Disconnect from power after use

14. TROUBLESHOOTING

Problem	Possible Cause	Solution
Does not turn on	No power	Check connection
Irregular flow	Clogged filter	Clean or replace
Air leakage from mask	Poor fit or deflated	Adjust and inflate with syringe (no needle)

15. TECHNICAL SPECIFICATIONS

Delivery rate: 10 LPM

Power consumption: 800W

Oxygen concentration: 92% max

Output pressure: 20–60 kPa

Deconcentrated oxygen: 9–16%

Voltage: 220V/50Hz – 110V/60Hz

Noise: <48 dBA

Weight: 20 kg

Output flow: 80–100 LPM

Dimensions: 34.4 × 30.6 × 56.5 cm

16. CONTRAINDICATIONS

PNEUMO[®]
PEAK
get your
own summit

Not recommended for individuals with:

Uncontrolled heart disease

Severe respiratory conditions

Pregnancy (unless medically indicated)

17. PROFESSIONAL USE AND LIABILITY

This equipment must be used under supervision of trained personnel. The manufacturer is not responsible for improper use.

18. WARRANTY AND TECHNICAL SERVICE

Warranty: 12 months

Covers manufacturing defects

Does not cover misuse or modifications

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Empresa: Pnemo Peak SAS
PNEUMOPEAK TECHNOLOGY S.A.S.